

Operator: John Doe  
 Process: Mapping  
 Session Start: 2010-05-12 15:28:01  
 Session Stop: 2010-05-19 09:41:01  
 Session Name: Environmental chamber mapping for web  
 Session Text:

Vessel: Environmental  
 Product: TSP loggers  
 Time Zone: W. Europe Daylight Time (UTC offset 02:00:00)



**Statistics**

Name	Statistics		
Description			
Type	Temperature		
Start Time	2010-05-12 15:28:01		
Stop Time	2010-05-19 09:39:01		
Global Min.	22.41 °C	2010-05-15 07:34:01	LA 02
		2010-05-15 07:35:01	LA 02
		2010-05-15 07:36:01	LA 02
		2010-05-15 07:37:01	LA 02
		2010-05-15 07:38:01	LA 02
		2010-05-15 07:39:01	LA 02
		2010-05-15 07:40:01	LA 02
		2010-05-15 07:41:01	LA 02
		2010-05-15 07:42:01	LA 02
		2010-05-15 07:43:01	LA 02
		2010-05-15 07:44:01	LA 02
		2010-05-15 07:45:01	LA 02
		2010-05-15 07:46:01	LA 02
		2010-05-15 07:47:01	LA 02
		2010-05-15 07:48:01	LA 02
		2010-05-15 07:49:01	LA 02
		2010-05-15 07:50:01	LA 02
		2010-05-15 07:51:01	LA 02
		2010-05-15 07:53:01	LA 02
		2010-05-15 07:55:01	LA 02
		2010-05-15 07:56:01	LA 02
		2010-05-15 07:57:01	LA 02
		2010-05-15 07:58:01	LA 02
		2010-05-15 07:59:01	LA 02
		2010-05-15 08:00:01	LA 02
		2010-05-15 08:02:01	LA 02
		2010-05-15 08:04:01	LA 02
		2010-05-15 08:05:01	LA 02
		2010-05-15 08:06:01	LA 02
		2010-05-15 08:07:01	LA 02
		2010-05-15 08:08:01	LA 02
		2010-05-15 08:09:01	LA 02
		2010-05-15 08:10:01	LA 02
		2010-05-15 08:11:01	LA 02
		2010-05-15 08:14:01	LA 02
Global Max.	24.59 °C	2010-05-17 20:18:01	LA 03
		2010-05-17 20:19:01	LA 03
		2010-05-17 20:20:01	LA 03
		2010-05-17 20:21:01	LA 03
		2010-05-17 20:22:01	LA 03
Global Average	22.89 °C		
Min average	22.85 °C	LA 02	
Max average	22.93 °C	LA 03	
Global Delta	2.18 °C		
Global Standard Deviation	0.35 °C		



Operator/Date: