

Operator: John Doe  
Process: Steam Sterilization  
Session Start: 14-08-2007 14:30:00  
Session Stop: 14-08-2007 17:59:00  
Session Name: Test  
Session Text:

Vessel: Autoclave 1  
Product: LVP  
Time Zone: WesternTime (UTC offset -04:00:00)



### Lethality Overview

Ref. Temperature: 121.11 °C      Z-Value: 10 °C  
Min. Temperature : 91.11 °C  
Start Date : 14-08-2007

Name:	Description:	Start Time:	End Time:	Max. Value:
LC 1 (F0) Lethality 1	L-Calculation	14:30:00	17:59:00	26.90
LC 2 (F0) Lethality 1	L-Calculation	14:30:00	17:59:00	27.75
LC 3 (F0) Lethality 1	L-Calculation	14:30:00	17:59:00	30.81
LC 4 (F0) Lethality 1	L-Calculation	14:30:00	17:59:00	27.88
LC 5 (F0) Lethality 1	L-Calculation	14:30:00	17:59:00	25.08
LC 6 (F0) Lethality 1	L-Calculation	14:30:00	17:59:00	22.68
LC 7 (F0) Lethality 1	L-Calculation	14:30:00	17:59:00	27.92
LC 8 (F0) Lethality 1	L-Calculation	14:30:00	17:59:00	26.05
LC 9 (F0) Lethality 1	L-Calculation	14:30:00	17:59:00	27.18
LC 10 (F0) Lethality 1	L-Calculation	14:30:00	17:59:00	26.75
LC 11 (F0) Lethality 1	L-Calculation	14:30:00	17:59:00	26.80

Operator/Date: