

Operator: John Doe
 Process: Heat Distrib
 Session Start: 2010-05-11 12:33:31 PM
 Session Stop: 2010-05-11 2:46:18 PM
 Session Name: Stability chamber test web 1
 Session Text:

Vessel: Stability Chamb
 Product: Loggers
 Time Zone: W. Europe Daylight Time (UTC offset 02:00:00)



Statistics

Name: Statistics on first run
 Description: Heat distribution
 Type: Temperature
 Start Time: 2010-05-11 12:40:06 PM
 Stop Time: 2010-05-11 1:09:24 PM
 Min average: 24.21 ° C LA 03
 Max average: 25.45 ° C LA 06

ID	Name	Description	Min.	Max.	Average	Delta	Standard Deviation
17235	LA 01	Temperature	24.47 ° C	26.45 ° C	25.23 ° C	1.98 ° C	0.52 ° C
21132-1	LC 02	Temperature	23.87 ° C	25.30 ° C	24.72 ° C	1.43 ° C	0.22 ° C
21132-2	LA 03	Temperature	23.76 ° C	25.31 ° C	24.21 ° C	1.55 ° C	0.37 ° C
21132-3	LA 04	Temperature	24.13 ° C	25.40 ° C	24.62 ° C	1.27 ° C	0.25 ° C
21132-4	LA 05	Temperature	23.93 ° C	25.33 ° C	24.36 ° C	1.40 ° C	0.29 ° C
17512	LA 06	Temperature	24.62 ° C	26.50 ° C	25.45 ° C	1.88 ° C	0.54 ° C
29257-2	LA 08	Temperature	24.52 ° C	25.55 ° C	25.03 ° C	1.03 ° C	0.33 ° C

SAMPLE

Operator/Date:

Operator: John Doe
 Process: Heat Distrib
 Session Start: 2010-05-11 12:33:31 PM
 Session Stop: 2010-05-11 2:46:18 PM
 Session Name: Stability chamber test web 1
 Session Text:

Vessel: Stability Chamb
 Product: Loggers
 Time Zone: W. Europe Daylight Time (UTC offset 02:00:00)



Statistics

Name: Statistics on second run
 Description: Heat distribution
 Type: Temperature
 Start Time: 2010-05-11 1:39:55 PM
 Stop Time: 2010-05-11 2:40:28 PM
 Min average: 23.02 ° C
 Max average: 24.92 ° C

LA 05
 LA 06

ID	Name	Description	Min.	Max.	Average	Delta	Standard Deviation
17235	LA 01	Temperature	24.05 ° C	26.21 ° C	24.89 ° C	2.16 ° C	0.60 ° C
21132-1	LC 02	Temperature	22.99 ° C	23.81 ° C	23.12 ° C	0.82 ° C	0.10 ° C
21132-2	LA 03	Temperature	23.01 ° C	23.85 ° C	23.16 ° C	0.84 ° C	0.09 ° C
21132-3	LA 04	Temperature	22.65 ° C	23.75 ° C	23.07 ° C	1.10 ° C	0.23 ° C
21132-4	LA 05	Temperature	22.66 ° C	23.81 ° C	23.02 ° C	1.15 ° C	0.27 ° C
17512	LA 06	Temperature	24.06 ° C	26.47 ° C	24.92 ° C	2.41 ° C	0.66 ° C
29257-2	LA 08	Temperature	23.93 ° C	26.57 ° C	24.77 ° C	2.64 ° C	0.72 ° C

SAMPLE

Operator/Date: