


Operator: John Doe  
 Process: P1 OP-SHORT-D  
 Session Start: 2007-10-17 10:02:28  
 Session Stop: 2007-10-17 10:45:41  
 Session Name: WD for Web application  
 Session Text: OQ Standard load 1 of 3

Vessel: WD Web  
 Product: EN 15883  
 Time Zone: W. Europe Daylight Time (UTC offset 02:00:00)



**Unit Configuration:**

Position:	Position Name:	Position Description:	Sample Rate:
1	1-1	Level 1 Right Top Corner	00:00:01
2	2-1	Level 1 Middle Centre	00:00:01
3	3-1	Level 1 Left Bottom Corner	00:00:01
4	4-1	Level 3 Left Top Corner	00:00:01
5	5-1	Level 3 Middle Centre	00:00:01
6	6-1	Level 3 Right Botom Corner	00:00:01
7	7-1	Level 5 Right Top Corner	00:00:01
8	8-1	Level 5 Middle Centre	00:00:01
9	9-1	Level 5 Left Bottom Corner	00:00:01
10	LAC	Automatic Control	00:00:01
11	HZ	Hot Zone	00:00:01
12	CZ	Cold Zone	00:00:01



Operator/Date: