

Operator: John Doe
 Process: P1 OP-SHORT-D
 Session Start: 2007-10-17 10:02:28
 Session Stop: 2007-10-17 10:45:41
 Session Name: WD for Web application
 Session Text: OQ Standard load 1 of 3

Vessel: WD Web
 Product: EN 15883
 Time Zone: W. Europe Daylight Time (UTC offset 02:00:00)



Statistics

Name	Statistics		
Description			
Type	Temperature		
Start Time	00:00:03		
Stop Time	00:42:56		
Global Min.	18,91 °C	2007-10-17 10:07:24	LA 6
Global Max.	93,71 °C	2007-10-17 10:31:55	LA 6
		2007-10-17 10:31:56	LA 6
		2007-10-17 10:31:57	LA 6
		2007-10-17 10:31:58	LA 6
		2007-10-17 10:31:59	LA 6
Global Average	56,88 °C		
Min average	55,62 °C	LA 5	
Max average	58,34 °C	LA 2	
Global Delta	74,80 °C		
Global Standard Deviation	19,47 °C		

ID	Name	Description	Min.	Max.	Average	Delta	Standard Deviation
21533-1	LA 1	Temperature	22,67 °C	93,61 °C	58,13 °C	70,94 °C	19,63 °C
21533-2	LA 2	Temperature	23,03 °C	93,62 °C	58,34 °C	70,59 °C	19,55 °C
21533-3	LA 3	Temperature	22,42 °C	93,59 °C	57,74 °C	71,17 °C	19,71 °C
21533-4	LA 4	Temperature	22,43 °C	93,58 °C	57,79 °C	71,15 °C	19,72 °C
21534-1	LA 5	Temperature	19,75 °C	93,67 °C	55,62 °C	73,92 °C	19,68 °C
21534-2	LA 6	Temperature	18,91 °C	93,71 °C	56,47 °C	74,80 °C	19,64 °C
21534-3	LA 7	Temperature	21,19 °C	93,65 °C	56,01 °C	72,46 °C	19,54 °C
21534-4	LA 8	Temperature	20,68 °C	93,66 °C	56,07 °C	72,98 °C	19,57 °C
21624-1	LA 9	Temperature	22,82 °C	93,61 °C	57,44 °C	70,79 °C	18,28 °C
21624-2	LA 10	Temperature	22,10 °C	93,63 °C	56,01 °C	71,53 °C	19,34 °C
21624-3	LA 11	Temperature	21,58 °C	93,65 °C	56,15 °C	72,07 °C	19,50 °C
21624-4	LA 12	Temperature	22,20 °C	93,66 °C	56,75 °C	71,46 °C	19,27 °C

Operator/Date: