RL1006 waterproof temperature transmitter

1000 Series

Part of the RL1000 range, the RL1006 waterproof temperature thermister is a sealed single plug in flexible thermister probe with 4pin lemo socket

Features

Dimension

- System complies with the European Radio Equipment Directive (RED) which covers EMC and radio specifications
- \checkmark A Battery Low alarm is displayed via PC software, warning of the requirement for battery change prior to power loss
- Software Alarm alert in the event of total transmission failure



Product code: RL1006-434.075 (other frequencies are available)

Radio transmitter functions	
Frequency options	A range of frequencies are available between 433-458MHz. Country specific regulations apply.
Radio power	10mW
Radio range	3km over open ground
Software required	W900 – Standard EMS Software Package W906 – Validated EMS Software Package *See EMS datasheet for further options
Hardware required	CR2 / CR3 – Controller SR2 – Smart Receiver REP – Repeater
Product codes	

Product codes	
88706	AA Lithium battery
Probe options	Please refer to the Hanwell Pro Sensor datasheet

Manufactured by Ellab Monitoring Solutions






Version 4- EC210078

Disclaimer: The information contained herein is believed to be reliable. Ellab Monitoring Solutions Ltd is not responsible for any incorrect or incomplete information on this datasheet and the information or product maybe changed without notice. Customers should obtain and verify the latest relevant information before placing orders for Hanwell products.

hanwell.com


Tel: +44 (0)1462 688070 | Email: sales.hanwell@ellab.com



UK CA

X

Œ



Instrumentation specification 167 x 40mm diameter 200 grams



