

# Wireless Environmental Monitoring Solutions



*Continuous monitoring with industry leading accuracy, reliability and radio range*

# Wireless Environmental Monitoring Solutions



Ellab's wireless environmental monitoring solutions provide highly accurate and reliable readings, while also providing the strongest signal coverage of any wireless system on the market today.

Our reputable wireless system is established globally within multiple industries including pharmaceutical, healthcare, food and industrial sectors. It can monitor a broad selection of applications such as fridges, freezers, chambers, incubators, labs, blood banks, manufacturing and research facilities, warehouse and storage by using the extended range of parameters.



Every customer's requirements are different, that's why we listen, ask questions and then offer our advice on the system that is right for your needs. We complete a Site Radio Survey of your sites; we commission your environmental monitoring systems and we offer you post-sales support for our products when you need it.






## Features

- A wide range of measurement options
- LCD display for a variety of transmitters and sensors
- Reliable, accurate and durable hardware suitable for long-term environmental monitoring
- Leading narrow band RF technology with transmission range up to 3,000 meters
- Long-life batteries that are user changeable
- 24/7 monitoring
- Automated alarms with various notification options
- Backup memory in transmitters and receivers
- Automated reporting
- FDA 21 CFR Part 11 compliant web-enabled EMS software





- 
**Accurate & reliable**  
 Quality assured hardware that is durable and easily maintained.
- 
**Secure**  
 Permission controls, data owned by customer.

- 
**Fully supported**  
 Service contract options available including hosted server option.
- 
**Assists with national and international regulatory compliance requirements**  
 FDA 21 CFR Part 11, GxP Guidelines, HACCP, EU GMP Annex 11 and ISO 22000 food standards.
- 
**Configurable to each user requirement**  
 Customized alarm notifications, data display, sensor & multi-site grouping.
- 
**Intuitive & audit ready**  
 Audit trails, scheduled reports, clear design for ease of use.
- 
**Validated software**  
 When required, we provide comprehensive IQ and OQ protocols and custom protocols for the Pharmaceutical sector. Trained engineers and distributors are available for smooth validation protocol execution.





# Hanwell Pro and IceSpy Wireless Solutions



Temperature



Differential Pressure



Humidity



Air Flow



CO<sub>2</sub>



Light & UV



Flood



Leak

## Hanwell Pro Solution

Pro is our superior wireless environmental monitoring system, that features high accuracy sensors, LCD display and unrivalled radio strength, which allows the system to overcome the most complex challenges.

Hanwell Pro wireless solutions provide advanced capabilities through an independent system capable of data analysis and reporting, which can be used to map and verify extensive ranges of temperature, relative humidity, light (LUX), UV, air flow, transport, damp, leak, wood, utilities and many more.

## Hanwell IceSpy Solution

IceSpy is our premium, reliable and cost effective wireless temperature monitoring system with highly competitive accuracy, cold-chain solutions and radio strength that overcomes most typical applications. For customers requiring validation systems, Hanwell IceSpy can also be fully validated to meet pharmaceutical, healthcare and food compliance requirements.


## Wireless technology you can trust

Nothing on the market can match the comprehensive technology that is available through the Hanwell Pro and IceSpy ranges.

The Hanwell Pro and IceSpy wireless monitoring systems embody decades of engineering expertise to ensure the best possible performance and guaranteed reliability. Combined with our unrivalled ability to plan, commission and execute robust systems and you have yourself a reliable partner.



# As your business needs grow - simply add more sensors



Option 1 Cloud Server  
Option 2 Internal Server

### Optional Alert Units

**1. ARB**

Alarm / RFID / Bluetooth

**3. Alarm Beacons**

Audio/visual alerts

**2. Autodialler**

Up to 9 call alerts

**4. SMS Module**

Out of hours alerts via text message

# Hanwell EMS

## - advanced software



Multisite/Group



Audit Trails



Historical Data



Email/SMS Alerts



Report Scheduling



RAG Alarms


Managing environmental data has never been so easy – whatever and whenever you need it – from wherever you are.

The highly comprehensive EMS software takes environmental monitoring data display to another level. Full data collection with interactive graphs, tables and plan views, allows users to slice historical environmental data in many ways for more advanced analysis.

EMS administrators can create user groups for tiered management of sensor groups, zones and/or multi-site management – making managing environmental data easier and more accessible than ever.

### Hanwell EMS Alarms App

The Hanwell EMS smartphone app is designed for users to access sensor alarm status, active alarms, list of sites, zones and sensor groups on-the-go.



Full overview



Plan view



Graphical view





**Whatever you need**



- Single or multi-user access for small, medium or large sites
- Share data with users and allocate permissions to various data groups
- Collect and display more data than any other system on the market
- Display quick overview of sites and/or sensor groups
- Access interactive graphs, tables and plan views for more detailed analysis

**Whenever you need it**



- Access real-time data, 24/7
- Immediate email, SMS and mobile alarm notifications to user-defined personnel
- Generate automated reports or alternatively access reports immediately and easy

**Wherever you are**



- Access critical information on your smartphone and local PC, via the intranet, network or internet.
- Locally installed software or cloud-based server hosting options are available
- View and interact with data and settings from anywhere in the world using your browser



**Calibration that provides peace of mind**

Our skilled calibration technicians provide flexible and cost-efficient maintenance of your investment and enable units to be calibrated either in-house or onsite.

We can test several calibration points and offer traceable and accredited ISO 17025 calibration and adjustment for temperature and humidity instruments.

Additionally, the Hanwell EMS software retains calibration information for each transmitter and automatically runs scheduled reports showing calibration due dates.

# Ellab



Ellab has been your validation and monitoring partner since 1949, offering wireless data loggers and thermocouple systems for thermal validation processes as well as wireless environmental monitoring solutions.

We serve both small and large companies within the life science and food industries and have solutions for almost all applications such as sterilization, freeze drying, environmental chamber testing, depyrogenation, warehouse mapping, pasteurization and many more.

Ellab develops unique and innovative solutions based on input from close interactions and dialogues with our customers. Our goal is to help our customers overcome challenges and increase their productivity by providing reliable and efficient solutions.

Ellab also offers complete turn-key or supplemental rental solutions, equipment qualification and validation services and our specialized training courses within Ellab Academy.

Ellab has a long tradition and dedication for delivering the highest performance and quality in our industry. Our user friendly and flexible validation and monitoring solutions are recognized and used by thousands of customers worldwide.

Ellab A/S is ISO 9001 & ATEX IEC 80079-34 certified. Our calibration laboratories in DK are accredited according to ISO/IEC 17025:2005 by DANAK under registration no. 520 and our German laboratory has a DAkkS ISO/IEC 17025:2005 accreditation. We also have a UKAS ISO/IEC 17025:2005 accreditation at our monitoring manufacturing site in UK.



Validation & Monitoring  
Solutions

## Ellab A/S

Trollesmindealle 25  
3400 Hilleroed  
Denmark  
+45 4452 0500

contact@ellab.com  
www.ellab.com

