

Operator: Validation Manager Vessel: Warehouse
 Process: Validation Product: Pharma
 Session Start: 22-11-2013 08:57:15 Time Zone: UTC offset 01:00:00
 Session Stop: 22-11-2013 08:57:25
 Session Name: Warehouse Validation
 Session Text:



Mean Kinetic Temperature (MKT)

Name: Mean Kinetic Temperature Report

Description:

Total Test Result: **Passed**

Input parameters

Activation energy (kJ/mol): 60
 Start time: 04-04-2007 22:00:00
 End time: 10-04-2007 16:00:00
 Time period (dd.hh:mm:ss): 00.23:00:00
 Number of time periods: 6

Success criteria

Minimum limit: 293,00 K 19,86 °C
 Maximum limit: 310,00 K 36,86 °C

Individual analysis

ID:	Sensor name:	Mean Kinetic Temperature:		Result:
90116	LC 14	293,51 K	20,37 °C	Passed
90127	LC 22	294,69 K	21,55 °C	Passed
90128	LC 41	293,56 K	20,42 °C	Passed
90129	LC 13	293,51 K	20,37 °C	Passed
90133	LC 8	293,53 K	20,39 °C	Passed
90134	LC 33	293,43 K	20,29 °C	Passed
90135	LC 42	294,28 K	21,14 °C	Passed
90136	LC 5	294,73 K	21,59 °C	Passed
90137	LC 25	293,20 K	20,06 °C	Passed
90138	LC 39	293,65 K	20,51 °C	Passed
90139	LC 46	294,89 K	21,75 °C	Passed
90140	LC 32	293,41 K	20,27 °C	Passed
90141	LC 6	293,59 K	20,45 °C	Passed
90142	LC 44	293,59 K	20,45 °C	Passed
90143	LC 2	294,69 K	21,55 °C	Passed
90144	LC 36	293,45 K	20,31 °C	Passed
90145	LC 4	294,55 K	21,41 °C	Passed
90147	LC 29	294,19 K	21,05 °C	Passed
90148	LC 23	293,78 K	20,64 °C	Passed
90149	LC 40	294,25 K	21,11 °C	Passed
90150	LC 38	294,10 K	20,96 °C	Passed
90151	LC 34	293,18 K	20,04 °C	Passed
90152	LC 11	294,31 K	21,17 °C	Passed
90153	LC 21	294,62 K	21,48 °C	Passed
90154	LC 28	293,14 K	20 °C	Passed
90155	LC 7	294,97 K	21,83 °C	Passed
90156	LC 27	293,95 K	20,81 °C	Passed
90157	LC 45	294,50 K	21,36 °C	Passed
90158	LC 48	293,50 K	20,36 °C	Passed
90159	LC 31	293,43 K	20,29 °C	Passed
90160	LC 17	294,83 K	21,69 °C	Passed

Operator: Validation Manager
Process: Validation
Session Start: 22-11-2013 08:57:15
Session Stop: 22-11-2013 08:57:25
Session Name: Warehouse Validation
Session Text:

Vessel: Warehouse
Product: Pharma
Time Zone: UTC offset 01:00:00



90161	LC 3	293,85 K	20,71 °C	Passed
90164	LC 26	293,30 K	20,16 °C	Passed
90165	LC 30	293,44 K	20,3 °C	Passed
90166	LC 35	294,27 K	21,13 °C	Passed
90167	LC 18	293,35 K	20,21 °C	Passed
90168	LC 20	293,78 K	20,64 °C	Passed
90169	LC 12	294,12 K	20,98 °C	Passed
90170	LC 19	293,39 K	20,25 °C	Passed
90171	LC 9	293,91 K	20,77 °C	Passed
90172	LC 1	294,05 K	20,91 °C	Passed
90173	LC 24	293,36 K	20,22 °C	Passed
90174	LC 37	293,36 K	20,22 °C	Passed
90176	LC 47	293,43 K	20,29 °C	Passed
90177	LC 43	293,26 K	20,12 °C	Passed
90178	LC 16	293,42 K	20,28 °C	Passed
90181	LC 10	293,72 K	20,58 °C	Passed
90182	LC 15	293,45 K	20,31 °C	Passed

SAMPLE