

Certification Issued Under the Authority of the Federal Communications Commission By:

TUV SUD BAPT UNLIMITED Octagon House, Concorde Way Segensworth North Fareham, PO15 5RL United Kingdom

Date of Grant: 11/28/2025 Application Dated: 11/26/2025

Ellab A/S Trollesmindealle 25 Hilleroed, 3400 Denmark

Attention: Gergely Sved , CEO

NOT TRANSFERABLE

EQUIPMENT AUTHORIZATION is hereby issued to the named GRANTEE, and is VALID ONLY for the equipment identified hereon for use under the Commission's Rules and Regulations listed below.

FCC IDENTIFIER: XUS-UXSKYDL1 Name of Grantee: Ellab A/S Equipment Class: Digital Transmission System Notes: The TrackSense Frigo UX SKY Loggers are part of a 2.4 GHz wireless measuring system with an Access Point connected to a PC with special software.

Table with 5 columns: Grant Notes, FCC Rule Parts, Frequency Range (MHZ), Output Watts, Frequency Tolerance, Emission Designator. Row 1: 15C, 2405.0 - 2475.0, 0.0025

Listed power is maximum rated. (measured is 1.9 mW peak conducted) RF Exposure exempt to routine evaluation at 0mm under 1.1307 (b)(3)(i)(A) 'Option A' after duty cycle correction of 0.35%, as evidenced in this filing.

